

**Dermatophagoides pteronyssinus
DERP1 Antibody, Rabbit PAb,
Antigen Affinity Purified**

Catalog Number: 69007-T16



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION

Immunogen:	Recombinant Dermatophagoides pteronyssinus DERP1 Protein (Catalog#69007-A08H)
Preparation	Produced in rabbits immunized with purified, recombinant Dermatophagoides pteronyssinus DERP1 (Catalog#69007-A08H; P08176; Met1-Leu320). DERP1 specific IgG was purified by Dermatophagoides pteronyssinus DERP1 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Dermatophagoides pteronyssinus DERP1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	DERP1

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Dermatophagoides pteronyssinus DERP1.
--------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.