

HGF/Hepatocyte Growth Factor Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 70063-R105

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Canine HGF/Hepatocyte Growth Factor Protein (Catalog#70063-DNAH)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Canine HGF/Hepatocyte Growth Factor (Catalog#70063-DNAH; NP_001002964.1; Met1-Ser730).
Ig Type:	Rabbit IgG
Clone ID:	105
Specificity:	Canine HGF/Hepatocyte Growth Factor
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Canine HGF/Hepatocyte Growth Factor.

Please Note: Optimal concentrations/dilutions should be determined by the end user.