

DLL1 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 80013-R006

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Rat DLL1 protein (Catalog#80013-R08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Rat DLL1 (rR DLL1; Catalog#80013-R08H; EDL99825.1; Met 1-Ser 534).
Ig Type:	Rabbit IgG
Clone ID:	006
Specificity:	Rat DLL1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Rat DLL1.

Please Note: Optimal concentrations/dilutions should be determined by the end user.