

B7-2/CD86 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 80025-R053

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Rat B7-2/CD86 Protein (Catalog#80025-R08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Rat B7-2/CD86 (rR B7-2/CD86; Catalog#80025-R08H; NP_064466.1; Met1-Ile250).
Ig Type:	Rabbit IgG
Clone ID:	053
Specificity:	Rat B7-2/CD86
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000
--------------	-----------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>