

IL-4R/CD124 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 80198-R002

EliteRmab® is a registered trademark of Sino Biological Inc.


GENERAL INFORMATION	
Immunogen:	Recombinant Rat IL-4R/CD124 Protein (Catalog#80198-R08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Rat IL-4R/CD124 (rR IL-4R/CD124; Catalog#80198-R08H; Q63257-1; Met1-Arg232).
Ig Type:	Rabbit IgG
Clone ID:	002
Specificity:	Rat IL-4R/CD124
Formulation:	0.2 µm filtered solution in PBS
APPLICATIONS	
Applications:	ELISA,FCM
RECOMMENDED CONCENTRATION	
Flow Cytometry	FCM: 1:25-1:100
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Rat IL-4R/CD124.

Please Note: Optimal concentrations/dilutions should be determined by the end user.



Catalog Number: 80198-R002

EliteRmab® is a registered trademark of Sino Biological Inc.



Flow cytometric analysis of rat IL4R(CD124) expression on SD rat splenocytes. Cells were stained with purified anti-rat IL4R(CD124), then a FITC-conjugated second step antibody. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.