

CD38 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 80229-T16

GENERAL INFORMATION

Immunogen:	Recombinant Rat CD38 Protein (Catalog#80229-R08H)
Preparation	Produced in rabbits immunized with purified, recombinant Rat CD38 (rR CD38; Catalog#80229-R08H; Q64244; Trp45-Val303). CD38 specific IgG was purified by Rat CD38 affinity chromatograph
Ig Type:	Rabbit IgG
Specificity:	Rat CD38
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
---------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000
-------	-----------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.