

CD6 / Cluster of Differentiation 6
Antibody, Rabbit PAb, Antigen
Affinity Purified

Catalog Number: 80311-T24



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION	
Immunogen:	Recombinant Rat CD6 / Cluster of Differentiation 6 protein (Catalog#80311-R08H)
Preparation	Produced in rabbits immunized with purified, recombinant Rat CD6 / Cluster of Differentiation 6 (Catalog#80311-R08H; Q812A4; Met1-Asn396). CD6 / Cluster of Differentiation 6 specific IgG was purified by Rat CD6 / Cluster of Differentiation 6 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Rat CD6 / Cluster of Differentiation 6
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA,IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:1000-1:4000
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Rat CD6 / Cluster of Differentiation 6.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CD6 / Cluster of Differentiation 6
Antibody, Rabbit PAb, Antigen
Affinity Purified


Catalog Number: 80311-T24




Sino Biological
Biological Solution Specialist



Immunohistochemical staining of rat CD6 in rat lymph node with rabbit polyclonal antibody (1:2000, formalin-fixed paraffin embedded sections).



Immunohistochemical staining of rat CD6 in rat spleen with rabbit polyclonal antibody (1:2000, formalin-fixed paraffin embedded sections).



Immunohistochemical staining of rat CD6 in rat thymus with rabbit polyclonal antibody (1:2000, formalin-fixed paraffin embedded sections).