

ART4 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 80383-T24

GENERAL INFORMATION

Immunogen:	80383-R08H
Preparation	Produced in rabbits immunized with purified, recombinant Rat ART4 (rR ART4; Catalog#80383-R08H; NP_001166980.1; Met1-Lys269). ART4 specific IgG was purified by Rat ART4 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Rat ART4
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	ARTC4

APPLICATIONS

Applications:	ELISA,IHC-P
---------------	-------------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:1000-1:4000
ELISA	ELISA: 1:5000-1:10000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

ART4 Antibody, Rabbit PAb, Antigen Affinity Purified

Catalog Number: 80383-T24



Sino Biological
Biological Solution Specialist



Immunohistochemical staining of rat ART4 in rat kidney with rabbit polyclonal antibody (1:2000, formalin-fixed paraffin embedded sections).