

Layilin / LAYN Antibody, Rabbit PAB



Sino Biological
Biological Solution Specialist

Catalog Number: 80431-RP01

GENERAL INFORMATION	
Immunogen:	Recombinant Rat Layilin / LAYN protein (Catalog#80431-R08H)
Preparation	Produced in rabbits immunized with purified, recombinant Rat Layilin / LAYN (Catalog#80431-R08H; NP_001178926.1; Met1-Glu224). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Rat Layilin / LAYN
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Rat Layilin / LAYN.

Please Note: Optimal concentrations/dilutions should be determined by the end user.