

CTHRC1 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 80712-R003

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Rat CTHRC1 Protein (Catalog#80712-R08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Rat CTHRC1 (rR CTHRC1; Catalog#80712-R08H; NP_001258229.1; Met1-Lys230).
Ig Type:	Rabbit IgG
Clone ID:	003
Specificity:	Rat CTHRC1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000
--------------	-----------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.