

Hemopexin / Hpx Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 81025-RP01

GENERAL INFORMATION	
Immunogen:	Recombinant Rat Hemopexin / Hpx Protein (Catalog#81025-R08H)
Preparation	Produced in rabbits immunized with purified, recombinant Rat Hemopexin / Hpx (rR Hemopexin / Hpx; Catalog#81025-R08H; P20059.3; Met1-Gln460). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Rat Hemopexin / Hpx
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Rat Hemopexin / Hpx.

Please Note: Optimal concentrations/dilutions should be determined by the end user.