

Influenza A Virus Hemagglutinin / HA Antibody (HRP), Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 86001-RM01-H

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant H1N1 HA protein
Preparation:	This antibody was obtained from a rabbit immunized with purified, recombinant Influenza A virus H1N1 Hemagglutinin extracellular domain and conjugated with horseradish-peroxidase (HRP).
Ig Type:	Rabbit IgG
Clone ID:	2
Specificity:	H1N1 (A/California/04/2009) HA H1N1 (A/California/07/2009) HA
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 0.1-1 µg/mL

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>