

IL18 / Interleukin-18 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 90011-MM02

GENERAL INFORMATION

Immunogen:	Recombinant Rhesus IL18 / Interleukin-18 Protein (Catalog#90011-C07E)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Rhesus IL18 / Interleukin-18 (Catalog #90011-C07E; NP_001028006.1; Ala 2-Asp 193). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG2a
Clone ID:	5G10E10
Specificity:	Rhesus IL18 / Interleukin-18
	No cross-reactivity with E.coli cell lysate in ELISA
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Rhesus IL18.
--------------	---

Please Note: Optimal concentrations/dilutions should be determined by the end user.