

CD32a / Fc gamma RIIA Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 90015-RP02

GENERAL INFORMATION	
Immunogen:	Recombinant Cynomolgus CD32a / Fc gamma RIIA protein (Catalog#90015-C08H)
Preparation	Produced in rabbits immunized with purified, recombinant Cynomolgus CD32a / Fc gamma RIIA (rC CD32a / Fc gamma RIIA; Catalog#90015-C08H; AAL92096.1; Met 1-Ile 211). CD32a / Fc gamma RIIA specific IgG was purified by Cynomolgus CD32a / Fc gamma RIIA affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Cynomolgus CD32a / Fc gamma RIIA
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Cynomolgus CD32A.

Please Note: Optimal concentrations/dilutions should be determined by the end user.