

HER3 / ErbB3 Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 90043-RP01

GENERAL INFORMATION	
Immunogen:	Recombinant Rhesus HER3 / ErbB3 protein (Catalog#90043-K08H)
Preparation	Produced in rabbits immunized with purified, recombinant Rhesus HER3 / ErbB3 (rh ERBB3; Catalog#90043-K08H; XP_005571238.1 / XP_001113953.2; Met1-Thr643). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Rhesus HER3 / ErbB3
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	ERBB3
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Rhesus HER3 / ErbB3.

Please Note: Optimal concentrations/dilutions should be determined by the end user.