

IFNA13 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 90108-T48

GENERAL INFORMATION	
Immunogen:	Recombinant Cynomolgus IFNA13 protein (Catalog#90108-C08H)
Preparation	Produced in rabbits immunized with purified, recombinant Cynomolgus IFNA13 (Catalog#90108-C08H; G7NFW4; Met1-Glu190). IFNA13 specific IgG was purified by Cynomolgus IFNA13 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Cynomolgus IFNA13
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	WB,ELISA
RECOMMENDED CONCENTRATION	
Western Blot	WB: 1:500-1:1000
ELISA	ELISA: 1:25000-1:50000 This antibody can be used at 1:25000-1:50000 with the appropriate secondary reagents to detect Cynomolgus IFNA13.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

IFNA13 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 90108-T48



Anti-IFNA13 rabbit polyclonal antibody at
1:500 dilution
Lane A: HL-60 Whole Cell Lysate

Lysates/proteins at 30 µg per lane.
Secondary
Goat Anti-Rabbit IgG H&L (Dylight800) at
1/10000 dilution.

Developed using the Odyssey technique.
Performed under reducing conditions.

Predicted band size: 22 kDa
Observed band size: 15 kDa
(We are unsure as to the identity of these
extra bands.)