

CD25/IL-2R alpha Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 90265-MM07

GENERAL INFORMATION

Immunogen:	Recombinant Cynomolgus CD25/IL-2R alpha Protein (Catalog#90265-C08H)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Cynomolgus CD25/IL-2R alpha (Catalog#90265-C08H; H6WS54; Met1-Arg213). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	07
Specificity:	Cynomolgus CD25/IL-2R alpha
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
---------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000
-------	----------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.