

ART1 / CD296 / ARTC1 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 90271-T24

GENERAL INFORMATION	
Immunogen:	Recombinant Cynomolgus ART1 / CD296 / ARTC1 Protein (Catalog#90271-C08H)
Preparation	Produced in rabbits immunized with purified, recombinant Cynomolgus ART1 / CD296 / ARTC1 (Catalog#90271-C08H; XP_005579037.1; Gln23-Asn294). ART1 / CD296 / ARTC1 specific IgG was purified by Cynomolgus ART1 / CD296 / ARTC1 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Cynomolgus ART1 / CD296 / ARTC1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA, IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:1000-1:4000
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Cynomolgus ART1 / CD296 / ARTC1.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

ART1 / CD296 / ARTC1 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 90271-T24



Immunohistochemical staining of cynomolgus ART1 in cynomolgus kidney with rabbit polyclonal antibody (1:2000, formalin-fixed paraffin embedded sections).