



Catalog Number: 90274-R101

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Rhesus CD4 Protein (Catalog#90274-C08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Rhesus CD4 (Catalog#90274-C08H; EHH20435.1; Met1-Trp390).
Ig Type:	Rabbit IgG
Clone ID:	101
Specificity:	Rhesus CD4
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA, IHC-P
----------------------	--------------

RECOMMENDED CONCENTRATION


ELISA	ELISA: 1:5000-1:10000
IHC-P	IHC-P: 1:50-1:200

Please Note: Optimal concentrations/dilutions should be determined by the end user.




Catalog Number: 90274-R101

EliteRmab® is a registered trademark of Sino Biological Inc.



Immunohistochemical staining of monkey CD4 in monkey spleen with rabbit polyclonal antibody at 1:100 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of monkey CD4 in monkey lymph node with rabbit polyclonal antibody at 1:100 dilution, formalin-fixed paraffin embedded sections.