

CD4 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 90274-T16

GENERAL INFORMATION

Immunogen:	Recombinant Rhesus CD4 Protein (Catalog#90274-C08H)
Preparation	Produced in rabbits immunized with purified, recombinant Rhesus CD4 (Catalog#90274-C08H; EHH20435.1; Met1-Trp390). CD4 specific IgG was purified by Rhesus CD4 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Rhesus CD4
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Rhesus CD4.
--------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.