

IL-1R1 / CD121a Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 90279-MM02

GENERAL INFORMATION	
Immunogen:	Recombinant Cynomolgus IL-1R1 / CD121a protein (Catalog#90279-C08H)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Cynomolgus IL-1R1 / CD121a (Catalog#90279-C08H; EHH61701.1; Met1-Thr332).
Ig Type:	Mouse IgG1
Clone ID:	02
Specificity:	Cynomolgus IL-1R1 / CD121a
	No cross-reactivity in ELISA with Cynomolgus Human IL1R2 Mouse IL1R1 Human IL1R1 Human cell lysate (293 cell line)
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Cynomolgus IL-1R1 / CD121a.

Please Note: Optimal concentrations/dilutions should be determined by the end user.