

c-Met / HGFR Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 90304-T16

GENERAL INFORMATION

Immunogen:	Recombinant Cynomolgus / Rhesus c-Met / HGFR protein (Catalog#90304-C08H)
Preparation	Produced in rabbits immunized with purified, recombinant Cynomolgus / Rhesus c-Met / HGFR (Catalog#90304-C08H; A0A2K5UM36-1 / NP_001162100.2; Met1-Thr932). c-Met / HGFR specific IgG was purified by Cynomolgus / Rhesus c-Met / HGFR affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Cynomolgus c-Met / HGFR
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	MET

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Cynomolgus / Rhesus c-Met / HGFR.
--------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>