

IL-23 (IL23A & IL12B Heterodimer) Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: CT012-RP01

GENERAL INFORMATION

Immunogen:	Recombinant Human IL-23 (IL23A & IL12B Heterodimer) protein (Catalog#CT012-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human IL-23 (IL23A & IL12B Heterodimer) (Catalog#CT012-H08H;NP_057668.1 & NP_002178.2; Met1-Pro189 & Met1-Ser328). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human IL-23 (IL23A & IL12B Heterodimer)
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	IL12B & IL23A

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human IL-23 (IL23A & IL12B Heterodimer).
--------------	---

Please Note: Optimal concentrations/dilutions should be determined by the end user.