

IL-23 (IL23A & IL12B Heterodimer) Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: CT028-R101

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Mouse IL-23 (IL23A & IL12B Heterodimer) Protein (Catalog#CT028-M08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Mouse IL-23 (IL23A & IL12B Heterodimer) (Catalog#CT028-M08H; Q9EQ14 & P43432; Met1-Ala196 & Met1-Ser335).
Ig Type:	Rabbit IgG
Clone ID:	101
Specificity:	Mouse IL-23 (IL23A & IL12B Heterodimer)
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse IL-23 (IL23A & IL12B Heterodimer).

Please Note: Optimal concentrations/dilutions should be determined by the end user.