

10×Digestion buffer 4 GMP-grade

Product description

This product is an optimized 10 × Digestion buffer for different restriction enzymes to ensure optimal(100%) digestion activity.The formulation of the product:500 mM Potassium Acetate,200 mM Tris-acetate,100 mM Magnesium Acetate,1mg/ml OsrHSA,pH7.9@25°C .

This product is produced in accordance with GMP process requirements and provided in a liquid form.

Specifications

Cat.No.	10668ES03 / 10668ES10 / 10668ES50
Size	5 T / 20 T / 50 T
Formulation	500 mM potassium acetate, 200 mM Tris-acetate, 100 mM magnesium acetate, 1 mg/ml OsrHSA, pH 7.9 at 25°C.

Components

Components No.	Name	10668ES03	10668ES10	10668ES50
10668	10×Digestion buffer 4 GMP-grade	1 mL	10 mL	50 mL

Storage

This product should be stored at -25~-15°C for 2 years.

Notes

1. Please operate with lab coats and disposable gloves, for your safety.