

10× Transcription Buffer 3

Product description

This product is optimized transcription buffer and efficient for in vitro synthesis of RNA. This product is mainly composed of buffer salts, reducing agents and magnesium ions required for enzyme proteins.

Specifications

Cat.No.	10674ES03 / 10674ES10 / 10674ES25/ 10674ES60 / 10674ES76
Size	1 mL / 10 mL / 25 mL / 100 mL / 500 mL

Components

Components No.	Name	10674ES03	10674ES10	10674ES25	10674ES60	10674ES76
10674	10× Transcription Buffer 3	1 mL	10 mL	25 mL	100 mL	500 mL

Storage

This product should be stored at -25~-15°C for one year.

Notes

1. This product is for research use only.
2. For your safety and health, please wear personal protective equipment (PPE), such as laboratory coats and disposable gloves, when operating with this product.