

Qubit Assay Tubes

Product description

Qubit Assay Tubes are disposable molecular labware designed for nucleic acid and protein quantification experiments. They are intended for use with mainstream Qubit fluorometers. These tubes are transparent, with the tube cap connected to the tube body in a single piece, featuring a snap cap design that provides an excellent user experience. Manufactured in a strictly controlled environment, these products are guaranteed to be sterile and free from DNase, RNase, and endotoxin contamination.

Specifications

Cat NO.	83509ES03 / 83509ES10
Size	500 tubes/box / 10 pcs/box

Features

1. Made from high-quality imported polypropylene material.
2. Compatible with Qubit fluorometers.
3. Thin and uniformly thick tube walls ensure good thermal conductivity.
4. The tube cap seals tightly when closed, preventing evaporation and contamination, and is easy to open.
5. Free of endotoxins, RNases, and DNases.

Storage

Ship and store at room temperature, valid life for 3 years.

Notes

1. For your safety and health, please wear a lab coat and disposable gloves during operation!
2. This product is for research use only!